## S<sub>N</sub>2 Reaction of an Alkyl Halide and an Aniline

#### CHEM HELP ASAP

experiment video: https://youtu.be/dFO2gI3Z1hQ

#### **Purpose**

The purpose of this experiment is to demonstrate an SN2 reaction between an alkyl halide (2-bromoacetophenone) and an aniline (para-toluidine).

### **Background**

SN2 reactions involve the attack of a strong nucleophile (often an anion) on the carbon of an alkyl halide. The nucleophile displaces the leaving group to break the carbon-halogen bond and form a new carbon-nucleophile bond. The reaction occurs in one step. If the starting halide is attached to a chiral carbon, then the stereochemical configuration is inverted during the reaction.

$$R \stackrel{\cdot \cdot \ominus}{\longrightarrow} R \stackrel{\cdot \cdot \ominus}{\longrightarrow} R$$

**Scheme 1.** Model S<sub>N</sub>2 reaction.

The  $S_N2$  reaction can be demonstrated on 2-bromoacetophenone (1) and *para*-toluidine (2) (Scheme 2). The reaction forms 2-[(4-methylphenyl)amino]-1-phenylethanone (3) after the included base, sodium bicarbonate, neutralizes the initial  $S_N2$  product.

**Scheme 2.** Today's reaction – an S<sub>N</sub>2 reaction

# Procedure – 2-[(4-methylphenyl)amino]-1-ethanone

In a 20 mL scintillation vial dissolve 2-bromoacetophenone (5.0 mmol) in 10 mL methanol with stirring. Add *p*-toluidine *trans*-stilbene (5.0 mmol) followed by sodium bicarbonate (6.0 mmol). All the reaction to stir in a shallow warm water bath. The temperature of the bath should not cause the reaction to boil. After 1 h check the reaction by TLC (30% EtOAc/70% hexane) to confirm consumption of either the halide or aniline. Cool the reaction mixture in an ice water bath. Filter the mixture with a Buchner funnel and a 125-mL side-arm flask. Be sure to seat the filter paper with methanol, not water. Rinse the filter cake with a minimal amount of cold methanol. Allow the isolated product to air dry, determine the product's mass, calculated a percent yield, record the melting range, and perform a TLC analysis of the product. Dispose of your TLC mobile phase in a waste container. Interpret the provided NMR spectrum. Record all your observations in your notebook.